




SARS-CoV-2 Complete Solutions

[Sample Collection • Extraction Kit • Detection Kit • Automations]




Sample Collection

Using Zybio Disposable Virus Sampling Tube




Nucleic acid Extraction and Purification

Using Pre-Pack Deep Well Plate Zybio Nucleic Acid Extraction Kit



SARS-CoV-2 Realtime PCR Detection

Using Zybio SARS-CoV-2 Detection Kit





Automation For Zybio Kit

EXM 3000: 9 Minutes Semi Automate Nucleic Acid Extractor



- Extraction up to 32 samples within 9 minutes
- Easy to operate with 5 programs that meet usual workload
- Safe and Reliable with negative pressure system, HEPA filtering and automatic UV sterilization

ZIP-96V: Qualitative Realtime PCR System

- User-friendly
easy to understand for standard usage
- Advance temperature control system
with 6 independent temperature control modules

